

Safety Scalpels

- Sterile
- Stainless steel, razor-sharp cutting edge
- Metric ruler on each scalpel
- Individually packaged

Reference	Size	Quantity
4-510	10	10/box
4-511	11	10/box
4-512	12	10/box
4-515	15	10/box
4-515C	15C Angled Blade	10/box
4-522	22	10/box



4-510



4-511



4-512



4-515



4-515C



4-522

Directions For Use

Gently push slider to right and completely forward in one simple motion. Blade locks into place with audible confirmation.
 To retract, push slider to right and back to return blade to original protected position.

