

Use Stoelting's Rodent Warmer X2 before, during and after surgical procedures to improve surgical outcome and overall longevity. The integrated touch screen allows for easy programming and controls 2 independent heating pads simultaneously. Available with a combination of mouse, rat or home-cage heating pads. **Rectal Thermal probe sold separately.**

Stoelting's Rodent Warmer X2 is ideal for use before, during or after surgical procedures. Warming the animal minimizes heat loss and improves surgical outcomes.



### Control two independent heating pads simultaneously

Minimize heat loss and improve surgical outcomes with Stoelting's Rodent Warmer X2.

#### Features:

- Integrated touch screen for easy programming
- Controls two independent heating pads simultaneously
- Rectal probe (sold separately) for body temperature feedback
- Available with a combination of mouse, rat or cage size heating pads

Heating pads are constructed of silicone rubber heaters.

Temperature Control Range	25-45°C
Temperature Resolution	0.1°C
Heater Blanket Connection	4-pin locking DIN
Heater Power	24VDC@3Amax
Power Requirements	120/240VAC (switchable), 50/60Hz, 75VA
Probe Input Connector	Phone jack
Probe Dimensions	1.66mm/0.62in tip diameter
Control Box Dimensions	14cm (l) x 13cm (w) x 4cm (h) - 5.5in x 5in x 1.6in
Control Box Weight	380 grams (13.4 ounces)
Mouse Heating Blanket Dimensions	7 x 7cm (2.75in x 2.75in)
Rat Heating Blanket Dimensions	15.25 x 15.25cm (6in x 6in)
Cage Heating Blanket Dimensions	16 x 38cm (6.25in x 15in)

Your Distributor is:  
 Braintree Scientific, Inc.  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)