


## Gas Anesthetizing Box



The 1/4 inch thick transparent plastic anesthetizing box allows for constant observation while anesthetizing. Helps protect against accidental overdosing.

For use with rats and mice. Has an easy to move machine tooled sliding top. Made of durable plastic, this box should last for years. Not for use with liquid organic solvents.

Available in the following sizes:	
Model	Description
AB - 1	(small box) 10"Lx4"Wx4"H
AB - 2	(large box) 10"Lx8"Wx8"H

Your Distributor is:  
Braintree Scientific, Inc.  
PO Box 850498, Braintree, MA 02185  
781-917-9526  
Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
Web: [www.braintreesci.com](http://www.braintreesci.com)