


Animal Rectal Probes The **New** isolated sensor (can be used with monitoring equipment where electrical interference from other electrodes may be present. Reinforced strain relief for applications where probe is in more frequent use. Smooth ball tip. Stainless steel shaft - 1" long, .059" dia. Tip diameter .125". 5 foot lead length. Time constant 0.8 secs.



Model	Dimension	
RET-2-ISO	Rectal probe for rats w/isolated sensor	
RET-3-ISO	Rectal probe for mice and neonatal rats w/isolated sensor	



The new RET-4 Rectal probe for very small animals (neonatal mice). Flexible polyurethane sheath with smooth rounded tip (see inset). 36" long (can be made longer). Outside diameter nominal 0.032" (max at tip 0.04"). Time constant 0.18 seconds. Max temp. 140°C. May be sterilized by Gas (ETO), Gluteraldehyde solution or autoclave.

RET-4	Rectal probe for very small animals (neonatal mice)	

Your Distributor is:
 Braintree Scientific, Inc.
 PO Box 850498, Braintree, MA 02185
 781-917-9526
 Email: Info@braintreesci.com
 Web: www.braintreesci.com