



Soft-Temp Veterinary Heating Pads

Conductive Warming For Veterinary Medicine

Soft-Temp reusable pads are the perfect solution for veterinary warming. They are heavy-duty and can be patched if needed. Available in 4 sizes, the pads fit table-tops, surgery tables and in kennel cages. Soft-Temp pads work best with Adroit's HTP-1500 pump, but adapters are available should you have a different pump*. Pads may be used directly against the skin or covered with pillow cases or light-weight fabric for added protection against damage.



V029 Medium



VKIT Patch Kit

Product Codes**

V016 Small 20x16 IN
V029 Medium 20x29 IN
V044 Large 20x44 IN
V057 X-Large 20x57 IN

VKIT Patch Kit
HTP-1500 Heat Pump

**Disposable pads and wraps also available, visit adroitmedical.com