


ATS-3


General Rodent Tattoo System

The most versatile lab animal tattoo identification system on the market !

**Lab Animal Tattoo Identification
Made Simple, Quick & Easy !**



- Heavy duty custom equipment.
- Humane mouse & rat restraints.
- Reduce animal stress (quick & safe).
- No anesthesia needed.
- Low cost (~5 cents/rodent).
- Detailed instructions (w/ DVD).
- After purchase support.
- Extended warranty and service programs.



ATS-3 Upgrade Options

Upgrade For Ultimate Flexibility & Convenience.

Neonate
Tattoo Identification

Option #1



Large Lab Animal
Tattoo Identification

Option #2



Battery Operation
Field Work Tattooing


Option #3



AIMS

Animal Tattoo Identification

Permanent, Easy, and Cost Effective



ATS-3 General Rodent Tattoo System has supplies and accessories for tattooing ~3000 mice/rats. The detailed training guide & instruction DVD provide step-by-step guidance. Use numbers, letters, lines, dots, and/or colors to permanently ID any color/strain mouse or rat. The most versatile, basic tattoo identification system for laboratory mice and rats.

ORDER: Cat #: ATS-3 General Rodent Tattoo System w/ supplies and accessories

Includes: Lab Animal Tattoo Identification Vol. 2. *Adult Mice and Rats*.

⇒ **Neonate Upgrade**

ORDER: Cat #: NEO-AP

Neonate Up-Grade Package

Includes: Lab Animal Tattoo Identification Vol. 1. *Neonate Rodents*

⇒ **Large Lab Animal Upgrade**

ORDER: Cat #: LLA-AP

Large Lab Animal Up-Grade Package

Includes: Lab Animal Tattoo Identification Vol. 3. *Large Laboratory Animals*

⇒ **Freedom Battery Upgrade**

ORDER: Cat #: BOP-AP

Freedom Battery Operation Up-Grade Package

Note: Requires ATS-3 Control Box Modification.

Upgrade
Options

All text and graphic are © 2009 by Animal Identification and Marking Systems, Inc.

All rights reserved. No text or graphics may be copied or used without the written permission from Animal Identification and Marking Systems, Inc.