

ANIMAL FEEDING NEEDLES

Flexible Teflon Needles - *New with Bite Protectors* (Disposable)


- Flexible Teflon tube is covered at the hub end with 304 stainless steel tubing – this protects tubing from damage due to biting during use
- Teflon tube is softer than other tubing types, very flexible, and translucent for easy confirmation of reagent delivery
- Plastic luer hubs and soft silicone tips
- Non-toxic and non-pyrogenic
- Packaged 5 needles per sterile Tyvek peel pouch - 20 pouches per box



Sterile Flexible Teflon Needles - *New With Bite Protectors* (Disposable)

Cat No.	Needle Style Material	Gauge X Length	Ball Diameter	Pack Size
BG-9928	Flexible Teflon PTFE Plastic w/ SS Bite Protector	18G X 1.5"	2.0mm	100
BG-9923		18G X 2"	2.0mm	100
BG-9924		18G X 3"	2.0mm	100
BG-9926		17G X 3.3"	2.4mm	100
BG-9925		15G X 3.1"	2.8mm	100
BG-9927		15G X 3.9"	2.8mm	100

**New With
Bite Protector**



Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com