



## Stainless Steel Brain Matrices

Your Distributor is:  
 Braintree Scientific, Inc.  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)



0.5mm Steel	1mm Steel	Size	Dimension A	Dimension B	Dimension C	Dimension D	Dimension E	Cavity Depth
BS-SS-505C	BS-SS-500C	Mouse coronal (40-75g)	3.18mm	11.1mm	8.73mm	19.1mm	12.2mm	7.4mm
BS-SS-505S	BS-SS-500S	Mouse sagittal (40-75g)	3.18mm	11.1mm	8.73mm	19.1mm	12.2mm	7.4mm
BS-SS-605C	BS-SS-600C	Small rat coronal (175-300g)	4.76mm	15.9mm	12.7mm	36.6mm	23.8mm	7.6mm
BS-SS-605S	BS-SS-600S	Small rat sagittal (175-300g)	4.76mm	15.9mm	12.7mm	36.6mm	23.8mm	7.6mm
BS-SS-705C	BS-SS-700C	Large rat coronal (300-600g)	4.76mm	19.8mm	14.7mm	36.6mm	24.7mm	10.9mm
BS-SS-705S	BS-SS-700S	Large rat sagittal (300-600g)	4.76mm	19.8mm	14.7mm	36.6mm	24.7mm	10.9mm