

LABNET INTERNATIONAL, INC.

MPS 1000™ Mini Plate Spinner

Finally, a compact and economical centrifuge for PCR plates

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com



- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges

Patent Pending



SPECIFICATIONS

Capacity:	2 PCR Plates
Rotor Speed:	2500 rpm
G-Force:	500 xg
Rotor:	Fixed Vertically
Dimensions (WxDxH):	7.5 x 8.2 x 7.2 in. 19 x 21 x 18.3 cm
Weight:	6 lbs / 3 kg
Electrical:	120V or 230V, 50-60 Hz

ORDERING INFORMATION

Catalog No.	Description
C1000*	Mini Plate Spinner centrifuge with two position rotor

* To order this product in 230V, please add (-230V) to the item number.

 **BRAINTREE
SCIENTIFIC, INC.**
LAB RESEARCH PRODUCTS

www.braintreesci.com / 781-917-9526