

Adam 54 5a_ bSuf EUS^W



The CB is a compact portable scale which is ideally suited for a variety of uses and applications.

It is small and easy to use. The bright green backlit LCD display shows clear results, with a simple to use keyboard and Grade 304 stainless steel removeable top pan for easy cleaning. Weigh units: grams, ounces and pounds, selectable.

The CB comes with batteries included and an optional mains adaptor can be purchased.

Model	Capacity	Increment	Pan Size
CB 501	1lb	0.0002lb	5.1"
CB 1001	2lb	0.0002lb	5.1"
CB 3000	6lb	0.002lb	5.1"

Optional: 6vDc 200mA Adaptor