

New! Checkit[®] Go – Measure dispensed volumes in just seconds

✓ Fast and Convenient

✓ Accurate

✓ Easy to Use

✓ Affordable

Engineered for Liquid Handling

Robots and Multi-Channel Pipettes

Breakthrough technology gets the job done in seconds!

Checkit Go – disposable cartridge instantly checks accuracy of liquid handlers and multi-channel pipettes

The Checkit Go is an easy-to-use disposable cartridge that allows you to check the accuracy of your liquid handling robot or multi-channel pipette in just seconds. No additional equipment is required.

Fast and Convenient – measure up to 8 channels in seconds

This surprisingly simple device measures the dispensed volumes of up to 8 channels simultaneously. It is the size of a well plate, so you can use it in your robot, or on a lab bench, or anywhere. You're done in under 10 seconds. That means no more trying to read a tiny droplet's weight on a high-precision balance, and repeat for 7 more channels.

Accurate – directly measures liquid volume

The Checkit Go directly measures the volume of the liquid dispensed, not an indirect measurement of the weight or the color of the liquid. In addition to being so quick and easy to use, the Checkit Go limits evaporation to negligible levels because the dispensed droplets are mostly encased in glass tubes within seconds. Measurements are guaranteed to be accurate to within 2% and typically better than 1%.

Your Distributor is:



www.braintreesci.com
781-917-9526