

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com

KEEPING ANIMALS HEALTHY, HAPPY & BEAUTIFUL

SUPER POCKET PRO® Compact Battery Operated Trimmer

CONVENIENTLY COMPACT

Wahl's **Super Pocket Pro®** battery-operated trimmer is the perfect size to store in a pocket when prepping for a show or doing quick touchups at home, ensuring your pet always looks its best when it matters most.

The Super Pocket Pro is manufactured with a durable chrome body and has been ergonomically designed to fit comfortably in your hand so you have complete control. This small, lightweight trimmer can be completely concealed in the palm of your hand for quick touch-ups on anxious animals.

The Super Pocket Pro kit includes a #40 blade and a #40T blade for precise, close-to-skin hair removal.



Features & Benefits:

- Quiet and easy to conceal in the palm of your hand
- BONUS: #40T Trimmer Blade and 6 Position Length Guide
- 30 Day Satisfaction Guarantee and 1 Year Limited Warranty

Ideal For:

- Trimming pet ears, faces, and paws
- Trimming cat hair

SUPER POCKET PRO® Compact Battery Operated Trimmer

1
YEAR
— LIMITED —
WARRANTY


30 **DAY**
RISK
FREE
100% SATISFACTION GUARANTEE

KIT CONTENTS



Battery Operated Compact Trimmer
#40 Trimmer Blade
BONUS: #40T Trimmer Blade
1 AA Battery
2 Attachment Guide Combs
6 Position Length Guide
Travel Pouch
Cleaning Brush, Blade Oil, and Instructions

SPECIFICATIONS



Model # 9961-2801
Color Black Chrome
UPC 0 43917 99618 9
SCC 0 43917 99618 6

Cord/Battery AA Battery
Speeds 1 Speed
SPM/RPM 6,000 SPM
Pack Style Blister
Case Pack 4
Manufactured China

Size (H x W x D) inches

Unit 4" x 1" x 1.5"
Pack 10" x 4.75" x 2.5"
Case 10" x 8" x 5.25"

Weight (oz. / lbs.)

Unit26 oz.
Pack39 lbs.
Case 1.87 lbs.