

CO₂ Chamber

for mice and rats

Door

- Place animals in chamber without losing CO₂.
- Provides for ultra quick animal induction.



Simple, Quick & Easy to Use!

The Cleanest CO₂ Induction Chamber on the Market!

Delivery Tube Within Chamber

- Fills chamber from bottom, decrease induction time, "more humane".

Quick Connect w/valve

- Closes when disconnected, keeps CO₂ in tube and chamber.

7' Tube

- Easily connects to CO₂ source.

Clear Acrylic Top

- Easy to see animals in chamber.



Solid Sides

- Prevents viewing of animals during the process.

Small Foot Print

- 9" diameter, easy to keep on bench.



Disposable Bag Liner

- Animal feces & urine remain in bag with animal.
- Chamber is always clean.

Cat #: COIC

CO₂ Chamber

Includes: See thru acrylic lid with door & quick connect w/ valve, 7' tubing, 50 disposable insert bags.



Cat #: COBL **Disposable Bag Liner (50 bags/pk)**

Use w/ CO₂ Chamber.

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com