

Product Information

The CurVet™ Rat Training Simulator

The Curvet rat training model is a purposely designed training aid to provide the skills required for humane and ethical handling and procedural competencies, without the need to use a live animal.

The Curvet rat model is British made using a material that does not include latex, but gives a unique sensitivity of a living specimen. This model is ideal for both veterinary and biomedical establishments and has several unique features such as:

- Realistic dermatological skin, suitable for:
 - ◇ Injections by most conventional routes
 - ◇ Micro-chipping identification
 - ◇ Handling and restraint (including the ability to scruff if required)
 - ◇ The deliberate absence of fur enables this model to be cleaned and sanitised with an appropriate chemical sterilant
- Head and spinal movement, including rear limbs
 - ◇ Suitable for oral administration
 - ◇ Lateral and ventral head movement
 - ◇ Hard palate with incisors
- Realistic removable tail with two lateral tail veins ideal for:
 - ◇ Blood sampling
 - ◇ Intravenous administration
 - ◇ Insertion of a flexible catheter (Angiocath)
 - ◇ Tail has a self sealing material hence long life span
 - ◇ Additional tails can be purchased separately
 - ◇ Artificial realistic non-toxic blood available
- Realistic features such as:
 - ◇ Eyes
 - ◇ Hard Palate with incisors
 - ◇ Pink extremities
 - ◇ Fore and rear limbs
 - ◇ Anus and rectum with limited aperture for thermoregulation monitoring i.e. clinical or digital thermometer
 - ◇ Adult rat size and weight approximately 300g

This model is NOT a toy, but a well thought out training aid to replace the use of a live animal to gain manual skill and competence, thus avoiding stress to both animal and handler.

Order Code RT001

The CurVet™ Rat Training Simulator



Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com