

ANIMAL FEEDING NEEDLES

Flexible Teflon Needles (Disposable)

- Teflon tube is softer than other tubing types, very flexible, and translucent for easy confirmation of reagent delivery
- Plastic luer hubs and soft silicone tips
- Non-toxic and non-pyrogenic
- Packaged 5 needles per sterile Tyvek peel pouch - 20 pouches per box



ANIMAL
FEEDING

Sterile Flexible Teflon Needles (Disposable)

Cat No.	Needle Style Material	Gauge X Length	Ball Diameter	Pack Size
DT-9931	Flexible Teflon PTFE Plastic	20G X 1.2"	1.9mm	100
DT-9920		20G X 1.5"	1.9mm	100
DT-9932		18G X 1"	2.0mm	100
DT-9933		18G X 1.2"	2.0mm	100
DT-9928		18G X 1.5"	2.0mm	100
DT-9923		18G X 2"	2.0mm	100
DT-9924		18G X 3"	2.0mm	100
DT-9922		17G X 2"	2.4mm	100
DT-9926		17G X 3.3"	2.4mm	100
DT-9925		15G X 3.1"	2.8mm	100
DT-9927		15G X 3.9"	2.8mm	100
DT-9929		14G X 4.3"	2.9mm	100
DT-9930		14G X 4.7"	2.9mm	100

Your Distributor is:
 Braintree Scientific, Inc.
 PO Box 850498, Braintree, MA 02185
 781-917-9526
 Email: Info@braintreesci.com
 Web: www.braintreesci.com