

“Blue Wax” Dissection Trays

Environmentally Safe Dissection Tray - **Completely Green**

“Blue Wax” can be heated to renew the surface and the dissection tray can be reused several times


“Blue Wax” is specially formulated to provide the proper consistency to allow dissection pins to stick firmly and be easily removed.

DTW1—Dissection Tray (9” x 6½” x 1”)



DTXL—X-Large Dissection Tray (12¼” x 8¼” x 1¼”)

“Blue Wax” dissection trays uses a completely environmentally safe wax product. Unlike the “Black Wax” **most commonly** used in other dissection trays. “Black Wax” dissection trays commonly list toxic chemicals and fire, health & reactivity hazards.



Use with AIMS
Procedure/Dissection Boards

Cat#: DTW1 AIMS Dissection Tray (9” x 6.5” x 1”)

Cat#: DTXL AIMS X-Large Dissection Tray (12¼” x 8¼” x 1¼”)