



Portable Animal Intensive Care Warmers

ThermoCare patented Portable Animal Intensive Care Warmer Units have been used successfully by laboratories, veterinarians, breeders and pet owners in aid of sick, recovering, or chronically ill animals, neonates and nursing mothers for over thirty years. Many original units, made in the late 1970's and early 1980's, are still in use today. The revolutionary heated lining of the water jacketed warmer, as well as the air within the warmer unit, are naturally circulated by convection currents. Safe, reliable, and versatile, ThermoCare's many solutions will meet the ever-changing needs of your diverse animal patients.

Benefits of the ThermoCare Portable ICU:

- **Safely Heat** post-operative patients and all cases of hypothermia without "hot spots" in a "safe heat" environment.
- **Humidify** the internal atmosphere; operating the humidifier from an oxygen line, air compressor or oxygen generator.
- **Nebulize** the internal atmosphere, without condensation; operated from an oxygen line or our optional air compressor.
- **Medicate** patients in lower respiratory stress with a nebulized mist of 1-3 micron size antibiotic particles.
- **Administer Oxygen** with a humidifier or nebulizer to moisten the atmosphere.
- **Birth, Brood or Whelp** hypothermic newborns with the "safe heat" of a water-jacketed warmer.
- **Incubate** eggs or newborns needing controlled heat therapy, oxygen and/or medicated nebulization after birth.


Size Class*	Model	Exterior Dimensions	Interior Dimensions
Small	#DW-1	27 1/2" L x 18" W x 17 1/2" H	21 1/2" L x 14" W x 14 3/4" H
Medium	#DW-2	33 1/8" L x 29 3/4" W x 26 1/2" H	30 3/4" L x 25" W x 24" H
Large	#DW-3	42 3/8" L x 29 3/4" W x 26 1/2" H	39 3/4" L x 25" W x 24" H
Small	#FW-1	27 1/2" L x 18" W x 8 3/4" H	21 1/2" L x 14" W x 6 1/4" H
Medium	#FW-2	33 1/4" L x 29 3/4" W x 14" H	30 3/4" L x 25" W x 12" H
Large	#FW-3	42 3/8" L x 29 3/4" W x 14" H	39 3/4" L x 25" W x 12" H



DW-1



FW-1



Accessories for use with Dome or Flat Cover

Oxygen Nebulizer, Medication Nebulizer, Humidifier, Oxygen Tubing, Hose Fittings and Connectors.

All Units Include: Digital Thermostat, Water Level Gauge, Spin/Venting Handport Covers, and a Humidigauge (temperature/humidity gauge).

Your Distributor is:

Braintree Scientific, Inc.

PO Box 850498, Braintree, MA 02185

781-917-9526

Email: Info@braintreesci.com

Web: www.braintreesci.com