

Elizabethan Collars

Elizabethan collars are designed to prevent animals from grooming and to protect surgical or treatment sites.

All collars are manufactured from a low-density polyethylene. All collars close with Velcro fasteners providing adjustable and secure fit for individual animals. Neck openings are bound with cotton jersey for comfort and light padding.



Elizabethan Collars

Catalog #	Size	Description	Weight Range	Interior Diameter		Exterior Diameter	
				in	cm	in	cm
MC E1		Mouse Elizabethan Collar		$\frac{3}{4}$	2.00	$2\frac{3}{4}$	7.00
RC E1		Rat Elizabethan Collar	up to 450g	$1\frac{1}{3}$	3.00	$4\frac{3}{4}$	12.00
RC E2		Rat Elizabethan Collar	450g+	2	5.00	$5\frac{1}{2}$	14.00
GP E1	Small	Guinea Pig Elizabethan Collar	up to 350g	$1\frac{5}{8}$	4.10	$5\frac{3}{8}$	13.72
GP E2	Medium	Guinea Pig Elizabethan Collar	350-600g	$1\frac{3}{4}$	4.45	6	15.24
GP E3	Large	Guinea Pig Elizabethan Collar	600g+	$1\frac{7}{8}$	4.76	$7\frac{1}{8}$	18.10
BC E0	X-Small	Rabbit Elizabethan Collar	2.7kg				
BC E1	Small	Rabbit Elizabethan Collar	2.7-3.6kg	$2\frac{1}{2}$	6.25	13	33.00
BC E2	Medium	Rabbit Elizabethan Collar	3.6-4.5kg	$2\frac{3}{4}$	7.00	$15\frac{1}{4}$	38.80
BC E3	Large	Rabbit Elizabethan Collar	4.5-5.4kg	3	7.70	16	40.70

P.O. Box 850498
Braintree MA 02185
www.braintresci.com

Phone (781) 917-9526
Fax (978) 244-8917
info@braintresci.com

