

EICKEMEYER® Dental Polisher / Drill Unit II

- Compact dental micro motor for straight and angled dental handpieces
- Powerful and durable, quiet operation
- Infinitely variable speed up to 35,000 rpm
- Forward and reverse movement
- Choice of hand or foot switch control
- Handpieces not included

Item No.

- | | |
|--------|---|
| 174028 | EICKEMEYER® Dental Polisher/Drill Unit II |
| 174026 | Low Speed Dental Handpiece – Contra Angle (RA) <ul style="list-style-type: none"> • For Dental Burr Set (174105) • Compatible with Dental Polisher & Drill Unit (174028) and EickSonic Plus (176880) • For polishing please order: Disposable Prophy Heads HP (174064) |
| 174027 | Low Speed Dental Handpiece – Straight (HP) <ul style="list-style-type: none"> • Ideal for rodent dentals • Use with long straight shank burrs, diamond discs and prophy heads • Compatible with Dental Polisher & Drill Unit (174028) and EickSonic Plus (176880) • For polishing please order: Disposable Prophy Heads HP (174064) |



174028



174026



174027

Dental Starter Kit for Small Animals

The most important equipment for EICKEMEYER® Dental Polisher/Drill Unit II (174028) and Dental Unit EickSonic Plus LED (176880) together in one set.

Item No.

- | | |
|----------|--|
| 174002 | Complete set, consisting of: |
| 174026 | Low Speed Dental Handpiece – Contra Angle (RA) <ul style="list-style-type: none"> • For Dental Burr Set (174105) • For polishing please order: Disposable Prophy Heads HP (174064) |
| 174027 | Low Speed Dental Handpiece – Straight (HP) <ul style="list-style-type: none"> • Ideal for rodent dentals • Use with long straight shank burrs, diamond discs and prophy heads • For polishing please order: Disposable Prophy Heads HP (174064) |
| 174105 | Dental Burr Set (RA) |
| 174064 | Disposable Polishing Heads (HP) <ul style="list-style-type: none"> • Single use • Max. speed 3,000 U/min • Compatible with Straight Angle (HP) Handpiece (174027) • Pack of 50 |
| 176781 | Prophylaxis Paste <ul style="list-style-type: none"> • Polishing paste • Content: 150 g |
| 174010 | Dental Handpiece Lubricant Oil <ul style="list-style-type: none"> • A specially formulated, highly effective general handpiece and air motor lubricant • Frequent application will assure hand piece and air motor performance and durability • Suitable for high & low speed handpieces • Content: 500 ml |
| 17401150 | Dental Handpiece Air Motor/Turbine Adaptor <ul style="list-style-type: none"> • For highspeed handpieces |
| 17401155 | Dental Handpiece Low Speed E-Type Adaptor <ul style="list-style-type: none"> • Suitable for 174010 |



174026



174027



174105



174064



176781



174010



17401150



17401155

Your Distributor is:
 Braintree Scientific, Inc.
 PO Box 850498, Braintree, MA 02185
 781-917-9526
 Email: Info@braintreesci.com
 Web: www.braintreesci.com