



**EZ-7000** - shown with E-Z Heat Controller, Heated Surgery Beds and Surgery Stations, EX-IC Exhaust Cube, and Chamber Warmer

## EZ-7000 CLASSIC SYSTEM

The E-Z Anesthesia Classic System is the redesigned version of the original E-Z Anesthesia System that has been used effectively for more than twenty years in hundreds of facilities.

- Accommodates animals ranging in size from neonate mice to rabbits
- Five outputs supply anesthesia to the induction chamber and up to four breathers
- Flow Meter configuration is more user-friendly
- Incorporates two flow meters:
  - Flow Meter 1) Controls flow into the vaporizer
  - Flow Meter 2) Controls flow into the induction chamber
- The second Flow Meter balances flow between the chambers and breathers
- Pushbutton oxygen flush purges induction chamber before removing animals
- Carry handle facilitates safe, easy transport of the unit

### Specifications:

Anesthetic: Isoflurane O - 5%  
Sevoflurane O - 8%

Carrier Gas: Oxygen or Air

Flow Rate: 0 - 4 LPM

Number of Outputs: 5

### EZ-7000 Base System Includes:

EZ-150C Classic Anesthesia Machine  
(1) Oxygen or Air Regulator  
(2) EZ-103A Microflex Breathers  
(1) EZ-178 Sure-Seal Mouse/Rat Induction Chamber  
(1) Case of ReFresh Charcoal Filters (8)

Also includes mouse/rat nosecones, nosecone plugs, plastic surgery bed, and all connecting tubes

### EZ-7000 Classic Base System



## OPTIONAL ACCESSORIES



### E-28000-TRAY

Optional Mayo Tray  
for accessories  
13.5" L x 9" W

*Equipment and  
accessories not  
included.*

### E-28000 Anesthesia Stand

- Accommodates each of our four anesthesia systems with minimal footprint
- Convenient, slide-on mounting
- Swivel top allows precise positioning
- Accommodates two 'E' size oxygen tanks
- Stable five-wheel base



### EZ-DOM Pole Mounted Dual O<sub>2</sub> Tank Manifold

### EX-IC EXHAUST CUBE

Optional EX-IC70 Exhaust Cube features an internal blower that captures waste gas from our vented chamber and up to four breathing devices.



Your Distributor is:  
Braintree Scientific, Inc.  
PO Box 850498, Braintree, MA 02185  
781-917-9526  
Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
Web: [www.braintreesci.com](http://www.braintreesci.com)