

## **EZ Clip™ Wound Closing Kit**

---

*For Rapid Suturing with Minimum Trauma*

**Ideal for use with lab animals or in veterinary surgery.** The EZ Clip™ kit consists of the 59023 EZ Clip™ Applier, case of 1000 Wound Clips (No. 59027) and the 59026 EZ Clip™ Remover.



### **EZ Clips™**

Racks of 9mm clips are individually packaged in heat sealed, peel-away envelopes for cleanliness and ease of handling and storage. EZ Clips™ are packaged 20 on a rack, with 10 racks per box (200 total), or 50 racks per case (1000 total). The clips are pyrogen free and non-sterile. However, the autoclavable plastic packaging allows for sterilization of the clips without removing them from the packaging.

### **EZ Clip™ Applier**

Stainless steel applying forceps hold one rack of 20 clips. Use as you would a staple gun to apply clips to wound. The clips are fed by rustless beryllium-copper springs.

### **EZ Clip™ Remover™**

Stainless steel EZ Clip™ Remover can be used to remove all sizes of wound clips.

Your Distributor is:  
Braintree Scientific, Inc.  
PO Box 850498, Braintree, MA 02185  
781-917-9526  
Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
Web: [www.braintreesci.com](http://www.braintreesci.com)