


*Grace
Boards*

*The New Standard for small animal surgery, procedure
& dissection platforms*

Procedure / Dissection Boards

Grace Boards are made of 1" thick High Density Polyethylene (HDPE) with excellent chemical resistance.




Magnetic Base

ORDER Cat#: **GPM2**




ORDER Cat#: **GPS2**




Magnetic Base

ORDER Cat#: **GPM1**



ORDER Cat#: **GPS1**



****Easy-to-Clean****
****Solid Work Surface****

**BRAINTREE
SCIENTIFIC, INC.**

Grace Procedure / Dissection Boards




Cat#	Description
GPM2	Grace Procedure/Dissection Board w/ LED, Magnifier, & Magnetic Base System Two Flexible LED lights - provides direct illumination with no shadows and no heat. Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3X, inset lens: 4.0X). (Allows continuous eye contact with object when moving into magnification area). Non-skid feet for added stability. Universal input power supply (100 – 240 VAC, 50-60 Hz). Size: 12" (L) x 9.5" (W).
GPS2	Grace Procedure/Dissection Board w/ LED, & Magnifier Two Flexible LED lights - provides direct illumination with no shadows and no heat. Flex arm, non-obstructive rimless magnifying glass (main lens: 2.3X, inset lens: 4.0X). (Allows continuous eye contact with object when moving into magnification area). Non-skid feet for added stability. Universal input power supply (100 – 240 VAC, 50-60 Hz). Size: 12" (L) x 9.5" (W).
GPM1	Grace Procedure/Dissection Board w/ Magnetic Base System Non-skid feet for added stability. Size: 12" (L) x 8" (W)
GPS1	Grace Procedure/Dissection Board Non-skid feet for added stability. 12" (L) x 8" (W).

Grace Boards are made of 1" thick High Density Polyethylene (HDPE) excellent chemical resistance.

Easy-to-Clean, Solid Work Surface

Accessories for Magnetic Base System

LED Light, Magnifier, Drape Holders, Limb Restraints



Large Magnetic Drape Holders

Set of 4

Cat#: **GDHL**



Small Magnetic Drape Holders

Set of 4

Cat#: **GDHS**



Magnetic Limb Restraints, Set of 4

Cat#: **GLR1**

Magnetic Ultra Bright LED Light

Universal input power supply
(100 – 240 VAC, 50-60 Hz)

Cat#: **GLED**

Magnetic Magnifier,
(main lens: 2.3X, inset lens:
4.0X). Cat#: **GMAG**

Blue Wax Dissection Tray
9" x 6.5"

Cat#: **GBDT**

 **BRAINTREE
SCIENTIFIC, INC.**

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com