

Illuminated Tail Injection Cone



“Illuminated Rodent Tail Vein Injection Cone”

The new 562-IRC features our “Bright Light” illumination system with rechargeable batteries and recharge unit, clear plastic tail insertion cone, formed white base with suction cup “non-slip” feet, housing batteries, on/off switch, lights, and an electrical plug for recharging the batteries (100-240 VAC capacity.)

Model Number

562-IRC

Cone Dimensions

5.75” long x 5” high x 4” opening

Animal Range

Mice: 15-30 grams
Rats: 500-600 grams

Your Distributor is:

Braintree Scientific, Inc.

PO Box 850498, Braintree, MA 02185

781-917-9526

Email: Info@braintreesci.com

Web: www.braintreesci.com