


## LabArt Mouse Restraint

Adjustable height restraint for laboratory mice. Provides optimal aid for injections and blood withdrawals from the tail vein. Restraint is made of transparent material to allow observations on the condition of the animal. The restraints are offered in 3 suitable sizes. Base plate autoclavable, Plexiglas restraint not autoclavable! Red restraint sleeve available to assist in calming of rodent.



### Features:

- Restraint consisting of Plexiglas with fixing ring and locking screw.
- Base plate with stainless steel guide rails and four slip-proof rubber feet.
- Base plate completely autoclavable.
- The restraint can be exchanged without any problems.

### Available in the following configurations:

Restraint Set - base plate and medium restraint (34 mm)

Restraint without base:

Restraint small (25 mm)

Restraint Medium (34 mm)

Restraint Large (54 mm)

Red Sleeve for restraint small

Red sleeve for restraint medium

Fully Autoclavable Restraint Set consisting of base plate and a medium restraint (34mm)

Autoclavable Restraint small (25 mm)

Autoclavable Restraint Medium (34 mm)



Your Distributor is:  
Braintree Scientific, Inc.  
PO Box 850498, Braintree, MA 02185  
781-917-9526  
Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
Web: [www.braintreesci.com](http://www.braintreesci.com)