

# Disposable Cages for Mice

ANIMAL WELFARE | UNBEATABLE VALUE | EMPLOYEE SATISFACTION



## BETTER CAGE, BETTER ANIMAL ENVIRONMENT

The Innocage® is designed to maximize the use of state-of-the-art airflow delivered by the Innorack®. The position of the supply and exhaust ports in the lid, in combination with the transversal intracage airflow, results a gentle, no-draft animal environment.

Third-party research shows our cages have 15x lower levels of ammonia than traditional systems and improve breeding ([click for the breeding study](#)).

Third-party research indicates that our cages have 15X lower levels of NH<sub>3</sub> than traditional caging and also feature a gentle animal environment.

*“Ammonia and Carbon Dioxide Concentrations in Disposable and Reusable Ventilated Mouse Cages,”* by Jerald Silverman, David W Bays, Sheldon F Copper, and Stephen P Baker.

Published in JAALAS Vol 47, March 2008.

READ ARTICLE HERE:

<http://www.disposablecages.com/disposable-cages-validation.php>



**81 in<sup>2</sup>**  
FLOOR SPACE

Innocage® IVC Mouse

## DESIGNED FOR EFFICIENT OPERATIONS

Innovive disposable mouse cages are designed for maximum transportation, handling and processing efficiency:

25 cages < 13 in

A stack of 25 Innovive mouse cages fits into a standard changing station.

50 < 5

A stack of 50 Innovive mouse cages is shorter than a stack of 5 traditional mouse cages.

1 pallet = 1200 changes

About **1200 Innovive cages change-outs** fit on a single pallet.

The dense stacking offered by Innovive disposable rodent caging **reduces the cost of transportation and storage**; it **decreases handling time** and **increases cage changeouts efficiency**.

## LIGHTWEIGHT and ERGONOMIC CAGES

At **8 oz**, the Innocage is the lightest cage on the market. Technicians from vivariums around the world have weighed in with their opinion: a lighter cage makes a big difference!

### Innocage® Mouse

11.6"L x 7"W x 5"H

81 in<sup>2</sup> floorspace

8.8 oz cage | 1.9 oz food tray

 **Made in USA**

**1 CAGE BOTTOM, OVER 150 OPTIONS**

The Innovive disposable rodent caging system features a universal cage bottom with multiple options for bedding, lids, ventilation, feeders and water.

This results in *over 150 ways to configure your rodent caging* to fit your specific research needs.

**THREE DELIVERY OPTIONS for Innovive Disposable Cage Bottoms:**



Innocage®, empty	Pre-bedded, corn cob	Pre-bedded, ALPHA-dri®
Innovive standard disposable cage bottom arrives ready for you to add the bedding of your choice.	Innovive standard disposable cage bottom arrives pre-filled with 1/8" corn cob and ready for use.	Our standard disposable cage bottom arrives pre-filled with 1/8" ALPHA-dri®. Perfect for GLP studies.

**FIVE LID OPTIONS:**



IVC Single Filter	IVC Dual Filter	IVC Containment	Static	Transport
<b>Single-Filter IVC</b> Lid designed for standard disposable IVC rodent housing applications.	<b>Dual-Filter IVC Lid</b> has filters that cover both ventilation ports.	<b>Dual HEPA Filter</b> IVC Lid designed for rodent containment.	<b>Static Housing Lid</b> is perfect for small studies, contingency caging & other scenarios when a rack isn't required.	No-water-bottle reservoir makes this <b>Transport Lid</b> ideal for moving mice short or long distances.

All Innovive disposable IVC mouse cage lids are compatible with the universal Innocage® bottom. The assembled Innocage® provides 81 in<sup>2</sup> of floor space and 5 inches of headroom.

All our caging components are delivered double-bagged and irradiated.

**TWO FEEDER OPTIONS:**



Stainless Steel	Standard
High-volume hybrid metal-plastic feeder with a stainless steel insert; compatible with all Innovive disposable cages.	High-volume hybrid metal-plastic feeder that features a zinc electroplated steel insert. Compatible with all our disposable cages.

*\*\*\*NOTE: For our zinc-plated hybrid feeder we use a process called bright zinc electroplating, a technology commonly used in residential water service pipes.*

**FOUR WATERING OPTIONS:**



Aquavive® Empty	Aquavive® Pre-filled	Autowatering	Water Pouches
Our 300 ml Aquavive® water bottle is delivered empty for filling at your facility.	Innovive standard Aquavive® water bottle can be delivered pre-filled with acidified water.	Innovive disposable rodent caging can be used with any standard autowatering solution.	Filling bags with water at your facility? No problem - the standard Innovive lid is a perfect fit.

Your Distributor is:  
 Braintree Scientific, Inc.  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)