

| Part #   | Animal Size     | Intravascular<br>Tip O.D. | Stylet | Total<br>Length | Connects<br>To<br>Needle | Material     | Includes |
|----------|-----------------|---------------------------|--------|-----------------|--------------------------|--------------|----------|
| MTV-1    | All Mice & Rats | 1 F                       | Yes    | 28cm            | 23ga                     | Polyurethane | Needle   |
| MTV-2    | All Mice & Rats | 1 F                       | No     | 61cm            | 23ga                     | Polyurethane | Needle   |
| MTV-2S   | All Mice & Rats | 1 F                       | Yes    | 61cm            | 23ga                     | Polyurethane | Needle   |
| MTV-3    | All Mice & Rats | 1 F                       | No     | 90cm            | 23ga                     | Polyurethane | Needle   |
| MTV-5    | All Mice & Rats | 1 F                       | No     | 162cm           | 23ga                     | Polyurethane | Needle   |
| MTV-BMS1 | All Mice & Rats | 1 F                       | Yes    | 14cm            | 23ga                     | Polyurethane | Needle   |
| MTV-11   | All Mice & Rats | 32ga                      | No     | 28cm            | 25ga                     | Polyurethane | Needle   |
| MTV-33   | All Mice & Rats | 1 F                       | Yes    | 28cm            | Male<br>Syringe          | Polyurethane | Stylet   |

Your Distributor is:  
 Braintree Scientific, Inc.  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)