

# ANIMAL FEEDING NEEDLES

## Animal Feeding Guide

ANIMAL FEEDING

Species	Weight Range	Gauge	Length	Ball Diameter	Disposable	Reusable
Mice	Up to 14 g	24G	1"	1.25 mm		Straight: 7900
	15-20 g	22G	1" - 1.5"	1.25 mm		Straight: 7901, 7920
	20-25 g	20G	1" - 3"	1.9 - 2.25 mm	Malleable: 9921, 9916 Flexible Teflon: 9931, 9920	Straight: 7921, 7902, 7903 Curved: 7910, 7911
	25-30 g	18G	1" - 2"	2.0 - 2.3 mm	Malleable: 9918 Flexible Teflon: 9932, 9933, 9928, 9923 w/ SS Bite Protector: 9928B, 9923B	Straight: 7922, 7904, 7905 Curved: 7912, 7913
	30-35 g	18G	2" - 3.1"	2.0 - 2.3 mm	Malleable: 9917 Flexible Teflon: 9923, 9924 w/ SS Bite Protector: 9923B, 9924B	Straight: 7905, 7906 Curved: 7913, 7914
Rats	50-75 g	20G	1" - 3"	1.9 - 2.25 mm	Malleable: 9921 Flexible Teflon: 9931, 9920	Straight: 7921, 7902 Curved: 7910
	75-120 g	18G	1" - 1.5"	2.0 - 2.25 mm	Flexible Teflon: 9932, 9933, 9928 w/ SS Bite Protector: 9928B	Straight: 7922, 7904 Curved: 7912
	100-200 g	18G	2" - 3.1"	2.0 - 2.3 mm	Malleable: 9918, 9917 Flexible Teflon: 9923, 9924, 9928 w/ SS Bite Protector: 9923B, 9924B, 9928B	Straight: 7905, 7906 Curved: 7913, 7914,
		17G	1" - 3.3"	2.4 mm	Flexible Teflon: 9922, 9926 w/ SS Bite Protector: 9922B, 9926B	
		16G	2" - 3"	3.0 mm		Straight: 7923, 7908 Curved: 7915, 7916
	150-300 g	16G - 15G	3" - 4"	2.8 - 3.0 mm	Malleable: 9915, 9919 Flexible Teflon: 9925, 9927 w/ SS Bite Protector: 9925B, 9927B	Straight: 7908 Curved: 7916, 7946
200-350 g	14G - 13G	3" - 5"	2.9 - 4.0 mm	Flexible Teflon: 9929, 9930	Straight: 7909, 7924 Curved: 7917	
Hamsters	60-200 g	18G	1.5" - 2"	2.25 - 2.3 mm	Malleable: 9918 Flexible Teflon: 9928 w/ SS Bite Protector: 9928B	Straight: 7905 Curved: 7913

Note: g=grams, kg=kilograms

Your Distributor is:  
**Braintree Scientific, Inc.**  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)

# ANIMAL FEEDING NEEDLES

## Animal Feeding Guide (Continued)

Species	Weight Range	Gauge	Length	Ball Diameter	Disposable	Reusable
Guinea Pigs	250-300 g	18G	1.5" - 2"	2.0 - 2.3 mm	Malleable: 9918 Flexible Teflon: 9928, 9923 w/ SS Bite Protector: 9928B, 9923B	Straight: 7904, 7905 Curved: 7912, 7913
		17G	2" - 3.3"	2.4 mm	Flexible Teflon: 9922, 9926 w/ SS Bite Protector: 9922B, 9926B	
	350-450 g	16G - 15G	3" - 4"	2.8 - 3.0 mm	Malleable: 9915, 9919 Flexible Teflon: 9925, 9927 w/ SS Bite Protector: 9925B, 9927B	Straight: 7908 Curved: 7916, 7946
	400-600 g	14G - 13G	3" - 5"	2.9 - 4.0 mm	Flexible Teflon: 9929, 9930	Straight: 7909, 7924 Curved: 7917
Rabbits	1-3 kg	16G - 15G	3" - 4"	2.8 - 3.0 mm	Malleable: 9915, 9919 Flexible Teflon: 9925, 9927 w/ SS Bite Protector: 9925B, 9927B	Straight: 7908 Curved: 7916, 7946
Primates	Up to 4.5 kg	12G	3.5"	6.25 mm		Curved: 7919

Note: g=grams, kg=kilograms

Your Distributor is:  
 Braintree Scientific, Inc.  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)