

## Non-Flexible Stainless Steel Needles (Reusable)

We also offer a wide variety of straight and curved animal feeding needles from 12 gauge to 24 gauge stainless tubing and from 1" to 4" in length. The precision assembly of stainless steel balls and type 304 stainless steel tubing ensures a quality needle suited for feeding and dosing without injury to sensitive tissues.

Features:

- Fixed shape from stainless steel tubing
- Stainless steel balls
- Can be cleaned and re-used as needed
- Packaged and sold in dozen quantity non-sterile



## Non-Flexible Stainless Steel Needles (Reusable)

Cat No.	Needle Style Material	Gauge X Length	Ball Diameter	Pack Size
N-020		22G X 1.5"	1.25mm	12
N-008	Non-Flexible - Straight 304 Stainless Steel	18G X 1"	2.25mm	12
N-009	Non-Flexible - Curved 304 Stainless Steel	16G X 4"	3mm	12



Your Distributor is:  
 Braintree Scientific, Inc.  
 PO Box 850498, Braintree, MA 02185  
 781-917-9526  
 Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
 Web: [www.braintreesci.com](http://www.braintreesci.com)