



## Pureline® OC4000

Scivena Scientific Anesthesia Innovations is proud to bring the revolutionary Pureline® OC4000 O2 concentrator to the biomedical research community for the first time. The revolutionary system provides on-demand oxygen to any existing anesthesia system. While similar machines have been in use in Europe for years, Scivena Scientific has perfected the technology, added safety features, and now is making this money-saving equipment available here. The Pureline® OC4000 O2 Concentrator is FDA approved, UL tested, and backed by the industry's best three year warranty.

## Product Information

---

- Ready for use with and Scivena Scientific Anesthesia System or your existing anesthesia equipment
- Extracts FDA-approved USP93% pure O2 from the air and delivers on demand- **up to 5 LPM**
- Uses the industry standard DISS male thread connection
- Includes 10' of high pressure O2 hose with female DISS end fittings
- FDA approved, UL tested
- Oxygen purity and flow safety alarms
- Quiet 40-dBA operation
- Made in the USA, state-of-the-art manufacturing with ISO quality and TUV SUV safety certifications
- Industry best 3 year warranty

Your Distributor is:  
Braintree Scientific, Inc.  
PO Box 850498, Braintree, MA 02185  
781-917-9526  
Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
Web: [www.braintreesci.com](http://www.braintreesci.com)