

Portable Digital Microscope / Camera

Equip your surgery platform with a portable digital microscope/camera

Ideal for small rodent surgery, demonstrations, dissections & much more !!



Capture

- a photo
- real time video
- time-lapsed video
- project image





General Features:

- Resolution: 640x480 pixels
- Interface: USB 2.0
- Magnification: 10~200X Continuous
- Sensor: 1/4" Color CMOS
- Video frame Rate: up to 30 frame-per-second
- Illumination: built-in 8 white LEDs
- CD Content: DinoCapture application software, User's manual
- Compatible Windows 2000/XP, Vista and most MAC OS
- The USB connection eliminates the need for interface cards or expensive software
- The software for the collection and processing of photos is included
- Software is easy and convenient to load and use
- Automatic or manual color balance and control
- Automatic or manual brightness control

Product Requirements:

- Windows 20Intel Pentium III 1GHz00/XP,
- Vista processor or faster
- 256 MB RAM
- CD-ROM drive
- 16-bit color display monitor or higher
- Available Hi-speed USB 2.0 Port



ORDER:

PDMC

Digital Microscope/Camera

Includes: Mounting Bracket (DMCB), DinoCapture software, USB 2.0 extension

DMCB

Bracket for Microscope/Camera (PDMC)

TSPCK

Thermo Surgery Platform/Camera Kit

Includes: Thermo Surgery Platform with gel-pads, Digital Microscope/Camera, Mounting Bracket (DMCB), DinoCapture Software, USB 2.0 Extension, Gel Pad Reactivation Unit.

SSPCK

Standard Surgery Platform/Camera Kit

Includes: Standard Surgery Platform, Digital Microscope/Camera, Mounting Bracket (DMCB), DinoCapture Software, USB 2.0 Extension.