

In today's laboratories consideration of all aspects of an animals' well-being and welfare contribute to the quality of scientific results. Lomir is pleased to present our Rodent Sling, a versatile restraining device suitable for a wide range of applications.

RODENT RESTRAINT SLING

Frame

Polypropylene base with supports made from stainless steel, one inch thick base material provides a robust support that does not move. Stainless supports and connectors held with convenient thumbs screws allow for height adjustments. The frame is easily cleaned and dismantles when not in use.



Restraints

Full body restraint is provided with the choice of either a sling suit or sling cover. Both restrainers allow for access to all limbs and the head and tail, all close with Velcro fasteners allowing for a quick and secure restraint and an adjustable fit for individual animals. Sling suits are sold in five sizes according to the weight range of the animals. The rodent sling cover is available in a single size for rats up to 500g. Custom sizes are available on request.



Slings can be used to securely restrain rodents during many procedures. Choice of restraint in many sizes and range of adjustment make the rodent sling a valuable asset to any laboratory setting.

Other small animal species may also be restrained using the same frame with suits and slings available for mice and guinea pigs.



Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com