

# Tail Vein Injection Made Easier!

AVAILABLE EXCLUSIVELY FROM BRAINTREE SCIENTIFIC



- ▶ Laboratory rodents are largely color blind, so a red tube looks like a comforting shelter
- ▶ Red tube is transparent, allowing investigator to observe animal inside for normal body position and breathing
- ▶ Improved comfort during acute experimental procedures without sacrificing safety or welfare

## Red Tube Tailveiner® Restrainer

**The first easy-to-use mouse and rat restrainer specifically designed for tail vein injections.**

Slotted tube allows a rodent to be quickly pulled in by the tail. The animal is then restrained by a plug which slides easily. The Tailveiner is the fastest restrainer yet for intravenous injections in small animals.

**Order your Braintree Scientific catalog today!**



 **BRAINTREE  
SCIENTIFIC, INC.**  
LAB RESEARCH PRODUCTS

Tel. 781-917-9526 • Fax: 978-244-8917 • Email: [info@BraintreeSci.com](mailto:info@BraintreeSci.com)

[braintreesci.com](http://braintreesci.com)

