

Part Description: Hallowell Rat & Mouse Intubation Pack



These items are for use in conjunction with the [Rodent WorkStand](#).

Each pack is designed to hold all the tools and supplies needed to perform quick, efficient and minimally traumatic orotracheal intubation of small rodents. The Rat Intubation Packs come with a specially designed and molded, [autoclavable intubation specula](#), [umbilical tape](#), [scissors](#), a [Lidocaine applicator](#), an [endotracheal tube guide wire](#), an [incisor loop](#), a [sleeve of cotton swabs](#), a [mirror](#), [five 14g catheters](#), and our CD Manual which contains videos, stills, and tutorials of how to perform the intubations. All of these items are packaged neatly and conveniently in a [13" x 8" x 2" polypropylene organizer box](#) with each compartment labeled for contents, reordering information and what the part is used for it can be stored most anywhere and always be ready for use.

[Mouse](#) or [Rat Intubation](#) Packs do not include the [otoscope](#) as there are three approved models to choose from. Lidocaine HCL Jelly is not included.

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com