

Schuco S5000 Nebulizer



Introducing the Schuco® S5000 Nebulizer ... Combining Schuco quality and reliability with an affordable price.

- Suppressed sound
- Enhanced airflow
- Easy to use and clean
- Our most cost-effective nebulizer

Complete with all that's needed to begin nebulizer therapy:

- Compact compressor
- Nebulizer cup
- Angled mouthpiece
- Oxygen tubing
- Five extra compressor filters
- 5-year limited warranty

Schuco[®] S5000

STANDARD FEATURES

- Quiet operation ... less than 60 dBA
- Compact and lightweight ... just 3-1/2 lbs.
- Easy to use and clean
- Convenient cord storage compartment
- Substantial compressor pressure ... 35-45 psig
- High compressor operating flow rate ... 7-11 lpm
- Therapeutic particle size ... 0.5 µm-5µm
- AC operation ... 120V 60 Hz
- 5-year limited warranty



61399



61400

ORDERING INFORMATION

S5000 Compressor, nebulizer cup, angled mouthpiece, oxygen tubing,
5 extra compressor filters

ACCESSORIES

S66-0126 Replacement Filters, 5/pkg. (20 pkgs./case)

64085 Adult Aerosol Mask with nebulizer and 7' clear tubing (50/case)

64095 Pediatric Aerosol Mask with nebulizer and 7' clear tubing (50/case)

61399 Replacement Nebulizer Kit: nebulizer cup, 7' clear tubing, mouthpiece, tee,
6" flex tube (50/case)

61400 Replacement Nebulizer Kit: nebulizer cup, 7' clear tubing, mouthpiece, tee (50/case)

88100 AC to DC Power Converter for use in car

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com