



The Solid Brass Mouse Restrainer is designed to comfortably restrain a mouse for tail vein injections or sampling.

**Features:**

- Size = 2" x 4"
- Weight = 2 lbs

**Benefits:**

- Able to warm for animal comfort and vasodilatation of the tail
- Solid Brass maintains heat for extended periods
- Weight provides a stable and secure placement on the lab bench
- Can be sterilized via autoclave, ETO, liquid, VHP
- Long lasting

Model# SAI-MBR

Your Distributor is:  
Braintree Scientific, Inc.  
PO Box 850498, Braintree, MA 02185  
781-917-9526  
Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)  
Web: [www.braintreesci.com](http://www.braintreesci.com)