

"Squeekums" is an incredible animal training alternative that allows students, lab techs, and handlers to learn how to handle a rodent with safety and confidence.

- Fully articulated and extremely realistic. Head, feet, and limbs move in a natural manner.
- The tail is detachable and disposable with the capability of IV access at the caudal vein site.
- Ear sets are replaceable for training ear tagging procedures.
- Comes with its own ABS hard case, IV accessories, artificial training blood, and instructions.

One of the most accurate representations of an animal! Future variations will have anatomically correct pharynx, trachea, stomach, peroral feeding, and orotracheal intubation access. "Squeekums" weighs the same as an average adult male rat.

