

Clinical Probes



RET-1 - Rectal probe for humans, neonatal, rabbits, and larger animals. Flexible, vinyl covered, soft tipped. Does not cause discomfort. Max. Temp. 90°C. Time constant 5.0 secs. 5 ft. lead. Isolated.

OT-1 - For fast reading oral use. Ball-tipped stainless steel shaft, stainless handle. 5 ft. lead. Max. Temp. 125°C. Time constant 0.8 secs. Not Isolated.


DSP-1 - Disposable sheaths for use with OT-1 and RET-1, per 1000.

SST-1 - Skin surface probe, to be taped on. 1/4" dia. 10 Kt. gold sensor disc, 5 ft. lead. Max. Temp. 90°C. Time constant 0.15 secs. Not Isolated.

SST-2 - as SST-1 with plastic handle. Isolated.

ESO-1 - Esophageal probe. Flexible vinyl covered 1/8" dia. Max. Temp. 90° C. Time constant 5.0 secs. 5 ft. lead. Isolated.

Animal Rectal Probes



RET-2 - Rectal probe for rats. For fast intermittent measurements. Smooth ball tip. Stainless steel shaft - 1" long, 0.59" dia., tip dia. .125". 5 ft. lead. Max Temp. 125°C. Time constant 0.8 secs. Not Isolated.

RET-3 - Rectal probe for mice, etc. As RET-2, with shaft 3/4" long. .028" dia. Tip diameter .065". Time constant 0.5 secs. Not Isolated.

Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com