


BETTER PRODUCTS. BETTER SCIENCE.



Screw Tops

### SCREW TOPS

- One-piece screw tops made of Ultem with stainless steel sipper tubes.
- Autoclavable to 300°F (150°C).
- Available with any of our 5/16" diameter tube.

Our most popular sizes are listed below:

Catalog #	Tube
PCST1.5*	1½" straight, open
PCST1.5TD*	1½" straight, ball point
PCST51BTD*	5"-1" bend, ball point
FSPCST1.5**	1½" straight, open
FSPCST1.5TD**	1½" straight, ball point
FSPCST41BTD**	4"-1" bend, ball point

\*For use with MST8, MST16 and MST32.

\*\*For use with PC9RH8.5, PC16RH8.5, FSPC-8HT & FSPC-16HT



ASSEMBLY  
 AVAILABLE

### STANDARD STOPPERS

Sampling of available sizes; call for sizes not listed.

- Available in rubber or autoclavable neoprene.

Composition	Size	Composition	Size
Rubber	5½	Rubber	9½
Rubber, flanged	5½	Rubber	10
Rubber	6	Rubber	10½
Rubber	6½	Neoprene	6½
Rubber	7	Neoprene	7
Rubber	7½	Neoprene	8½
Rubber	8	Neoprene	10
Rubber	8½	Neoprene	10½

### DOUBLE SEAL STOPPERS

- Fully autoclavable to 300°F (150°C).
- Highly recommended for use with glass and Macrolon bottles.
- Equivalent Standard Stopper size is No. 6½.
- Prices are per stopper.



Double Seal Stoppers

Catalog #	Tube Size
S-100	5/16"
S-101	5/8"