

SPECIFICATION

Accuracy:	± 0.5°C
Resolution:	±0.1°C
Range:	-100 to +400°C
Readout:	Dual display of Set(green) and Run(red) temperatures
Max AC Power Output with resistive load:	TCAT-2DF and TCAT-2AC: 700 watts (Does not apply to TCAT-2LV)
DC Output Rating:	TCAT-2DF and TCAT-2LV: 12V DC, 2 Amps Max. Custom output voltages available upon request (Does not apply to TCAT-2AC)
AC ground leakage current:	less than 1 millamp at 110 V AC
Dual Setpoint Protection Feature:	TCAT-2DF and TCAT-2LV only additional heating pad feedback sensor prevents temperature exceeding 40°C..
Dimensions:	7 1/4" W x 4 3/8" H x 10" D
Weight :	Less than 5 lbs.
Power Supply:	110V or 220 V AC