

# T-FLOLOC™ - 2% Tetra-EDTA Flush and Lock Solution

*the PREMIUM catheter flushing and locking solution offering both patency and infection control*

## *Provides Broad-Spectrum Activity*

effective against gram positive & gram-negative bacteria and fungi associated with catheter-related infections

## *Prevents and Eradicates Biofilm*

eliminates protected microbial colonization of catheters effectively reducing the risk of catheter-related infection

## *Effective Anticoagulant*

promotes catheter patency; avoids systemic anticoagulation associated with heparin

## *Favorable Safety Profile*

non-toxic; no topical or systemic effects are known

## *No Development of Resistance*

does not induce development of drug-resistant bacterial strains



## T-FLOLOC™ ORDERING INFORMATION

Catalog Number	Description
TFL-1202-01	<u>1 Box</u> of 10 x 2mL 2% Tetra-EDTA Solution in 12mL Syringe
TFL-1202-05	<u>5 Box</u> of 10 x 2mL 2% Tetra-EDTA Solution in 12mL Syringe
TFL-1202-10	<u>10 Box</u> of 10 x 2mL 2% Tetra-EDTA Solution in 12mL Syringe

*NOTE: each syringe is sterile packaged, suitable for sterile field use*