

SMALL ANIMAL HANDLING

Gloves are suitable for the safe handling of small laboratory animals including

rodents. T/FLEX® PLUS gloves are a 13 gauge seamless

knit shell; material is Dyneema® engineered fibers plated with Lycra®

- Spandex®, Thermax® and the antimicrobial treatment

Alphasan® combine to provide a comfortable glove featuring

the highest cut resistance. Gloves conform to international

safety standards for cut resistance.



SAFE, FLEXIBLE AND COMFORTABLE

Available in five sizes: x-small through x-large

Cat #	Description
TEL-13G6	T Flex Plus handling gloves, x-small
TFL-13G7	T Flex Plus handling gloves, small
TEL-13G8	T Flex Plus handling gloves, medium
TFL-13G9	T Flex Plus handling gloves, large
TEL13-G10	T Flex Plus handling gloves, x-large

HANDLING MADE EASY



Your Distributor is:
 Braintree Scientific, Inc.
 PO Box 850498, Braintree, MA 02185
 781-917-9526
 Email: Info@braintreesci.com
 Web: www.braintreesci.com