

T/Pump® Localized Therapy System



Your Distributor is

**BRAINTREE
SCIENTIFIC, INC.**

Website: www.braintreesci.com • Email: info@braintreesci.com
Phone: 781-348-0768 • Fax: 781-843-7932

Standard Features

- 8.5' power cord
- Hospital grade plug
- 93 oz (2750 mL) reservoir
- Large fill port (1-7/8" diameter)
- Wet or dry storage

Optional Accessories

- TPS2 Temperature therapy stand with storage basket and cord hooks
- TPIV IV pole mount plate

Optional Test Accessories

- TFCI Test Thermometer
- TPT9 Flow Meter
- Test Probe Assembly

Connector Options

- Klik-Tite Connectors are leakproof and multiple pads can be connected to one unit for multiple therapy sites
- 600 Series Mul-T-Pads offer quick disconnect Colder Style Connectors for increased adaptability with TP650C TP700C models

Optional Mul-T-Pads

Single Patient Use with Klik-Tite Connectors

- TP22B 15" x 22" (38 x 56 cm)
Pliable all-polymer (set of 10)
- TP22C 15" x 22" (38 x 56 cm)
Pliable polymer/non woven fabric on one side

Lightweight, Single Patient Use with Klik-Tite connectors

- TP3E 3" x 23" (8 x 59 cm)
Pliable polymer/nonwoven fabric on one side
- TP12E 13" x 18" (33 x 46 cm)
Pliable polymer/nonwoven fabric on one side

Specifications

Model Number	
T/Pump Professional	TP700, TP700C*, TPP722**
T/Pump Classic	TP650, TP650C*, TPP622**
Height	11.5" (29 cm)
Width	8" (20 cm)
Depth	8" (20 cm)
Weight	
Full (filled to heating water line)	9.3 lbs (4.2 kg)
Empty	6.7 lbs (3 kg)
Temperature Range	10°C-42°C (50°F-107.°F)
Accuracy	±2°F (1°C)
Fluid System	9 gph flow rate min.

Electronics

Volts	120 VAC, 3.1 amps, 60 Hz, 20µA typical leakage current
Safety	UL approved to UL60601-1 and ASTM F2196
	EMC approved to EN60601-1-2, Class B

Warranty

1 year

Stryker reserves the right to change specifications without notice.

*TP700 and TP650 have Klik-Tite Connectors. TP700C and TP650C have Colder style Connectors.

**TPP722 and TPP622 come with the TP22C Mul-T-Pad.

- TP22E 15" x 22" (38 x 56 cm)
Pliable polymer/nonwoven fabric on one side

- TP26E 18" x 26" (46 x 66 cm)
Pliable polymer/non woven fabric on one side

Reusable

- TP22G 15" x 22" (38 x 56 cm)
Heavy duty pliable polymer

600 Series, Single Patient Use With Colder Style Connectors

- TP612E 13" x 18" (33 x 46 cm)
Polymer/nonwoven fabric on both sides
- TP622E 15" x 22" (38 x 56 cm)
Polymer/nonwoven fabric on both sides
- TP626E 18" x 26" (46 x 66 cm)
Polymer/nonwoven fabric on both sides