



## Trocars

Skin tunneling needles (trocars) offer a quick and easy method for routing catheters and electrodes subcutaneously. Our trocars are manufactured from thick walled medical grade stainless steel and have a semi-blunt 45 degree beveled tip or a stylet. The 7G trocar is perfect for most large animal catheters while the 10G is suited for smaller animals and catheters. Call for custom lengths.

TRO-7-6	Semi-Blunt Tip, 7 gauge x 6"/15cm
TRO-7-12	Semi-Blunt Tip, 7 gauge x 12"/30cm
TRO-10-8	Semi-Blunt Tip, 10 gauge x 8"/20cm

Trocar with stylet pictured above

## Trocars with Stylet

TRO-14MS - Mouse Trocar with Stylet

TRO-7B-12 - Semi-Blunt Tip, 7 gauge x 12 in. with Stylet

TRO-7B-6 Semi-Blunt Tip, 7 gauge x 6 in. with Stylet

Your Distributor is:

Braintree Scientific, Inc.

PO Box 850498, Braintree, MA 02185

781-917-9526

Email: [Info@braintreesci.com](mailto:Info@braintreesci.com)

Web: [www.braintreesci.com](http://www.braintreesci.com)