



Trocars

Skin tunneling needles (trocars) offer a quick and easy method for routing catheters and electrodes subcutaneously. Our trocars are manufactured from thick walled medical grade stainless steel and have a semi-blunt 45 degree beveled tip or a stylet. The 7G trocar is perfect for most large animal catheters while the 10G is suited for smaller animals and catheters. Call for custom lengths.

TRO-7-6	Semi-Blunt Tip, 7 gauge x 6"/15cm
TRO-7-12	Semi-Blunt Tip, 7 gauge x 12"/30cm
TRO-10-8	Semi-Blunt Tip, 10 gauge x 8"/20cm

Trocar with stylet pictured above

Trocars with Stylet

TRO-14MS - Mouse Trocar with Stylet
TRO-7B-12 - Semi-Blunt Tip, 7 gauge x 12 in. with Stylet
TRO-7B-6 Semi-Blunt Tip, 7 gauge x 6 in. with Stylet

P.O. Box 850498
Braintree MA 02185
www.braintresci.com

Phone (781) 917-9526
Fax (978) 244-8917
info@braintresci.com