



For many years, Radio Frequency Identification has provided laboratories with a secure method to identify animals. Unified Information Devices (UID) provides implantable transponders, readers, and multiple tag programming that is compatible with any LIMS system. UID also creates customized software solutions based on your individual laboratory processes in order to computerize your Excel sheets. UID transponders and readers integrate flawlessly with Provantis®, Pristima™, PDS, Topaz, and many others for easy and secure data collection. UID also has a complete line of Trovan and ISO Standard compatible transponders available.

IMPLANTABLE RFID TRANSPONDERS FOR ANIMALS


All UID transponders for animal research are compatible with the URH1 handheld reader. Each transponder is supplied in a triple ground cannula for easy insertion. Perfect for all animals from mice to horses. All are individually packaged and sterilized. For non-Sterile or bulk purchase, contact UID.

- Easy to insert
- Provided with a unique identification number.
Never to be repeated
- Can be programed using the UPS-60 Provisioning Station
for your specific 15 digit Identification



IMPLANTABLE RFID CHIPS FOR MICE

- Smallest Transponder on the market
- 1.4mm x 8.5mm in a sharp cannula
- Sterile and individually packaged



Product #	Description	Size (mm)	Qty
UC-1485	Mouse Glass (w/Cannula)	1.4 x 8.5	24
UC-1485Q	Mouse Glass (w/Cannula)	1.4 x 8.5	200
UC-1485R	Mouse Glass (w/Cannula)	1.4 x 8.5	1000

IMPLANTABLE RFID CHIPS FOR RATS

- The most widely used size for rodents and other animals
- 2.1mm x 12mm in a sharp cannula
- Sterile and individually packaged



Product #	Description	Size (mm)	Qty
UC-2112	Rat Glass (w/Cannula)	2.12 x 11.5	24
UC-2112Q	Rat Glass (w/Cannula)	2.12 x 11.5	200
UC-2112R	Rat Glass (w/Cannula)	2.12 x 11.5	1000

Your Distributor is:
Braintree Scientific, Inc.
 PO Box 850498, Braintree, MA 02185
 781-917-9526
 Email: Info@braintreesci.com
 Web: www.braintreesci.com

TRANSPONDER READER



UID readers are designed to work with the secure UID Transponders. Each transponder has an encryption that only allows for communication with a UID compatible reader. All UID Readers are fully compatible with any LIMS or other software that allows the input of Identification Data. Stand ISO readers are also available.

Hand Held Reader

- USB or Bluetooth connectivity
- Battery powered
- Available Protective case

Part Number: URH1



Available protective cases
(Available in orange or black)

CHIP PROVISIONING STATION

The most advanced multi-transponder programmer available. Our new provisioning station allows you to program up to 6 transponders at one time. Program the exact number or letter sequence specific to your needs.

Patent Pending

- Allows the programming of up to 6 UID injectable transponders at one time
- Code up to 15 alpha-numeric characters of your choice prior to implantation
- Fast and secure transponder programming
- See what you programmed on the reader when you scan the animal!

Part Number: UPS-60



Your Distributor is:
Braintree Scientific, Inc.
PO Box 850498, Braintree, MA 02185
781-917-9526
Email: Info@braintreesci.com
Web: www.braintreesci.com