

ANIMAL FEEDING NEEDLES

Flexible Teflon Needles (Disposable)

- Teflon tube is softer than other tubing types, very flexible, and translucent for easy confirmation of reagent delivery
- Plastic luer hubs and soft silicone tips
- Non-toxic and non-pyrogenic
- Packaged 5 needles per sterile pouch
20 pouches per box



Sterile Flexible Teflon Needles (Disposable)

Cat No.	Needle Style Material	Gauge X Length	Ball Diameter	Pack Size
DT-9920S	Flexible Teflon PTFE Plastic	20G X 1.2"	1.9mm	100
DT-9920		20G X 1.5"	1.9mm	100
DT-9928		18G X 1.5"	2.0mm	100
DT-7202		18G X 2"	2.0mm	100
DT-5202S		18G X 1.2"	2.0mm	100
DT-9926		17G X 3.3"	2.4mm	100
DT-9925		15G X 3.1"	2.8mm	100