


The patent pending STEREODRAFT is the first heated surgical table with bi-level vacuum designed specifically for standard stereotaxic frames using gas anesthesia. The bi-level vacuum, independent studies show the Stereodraft exceeds OSHA's requirements, will remove escaped gases from the atmosphere ensuring operator safety while maintaining proper animal temperature. Made of durable plastic that can withstand temperatures up to 140 degrees.

www.STEREODRAFT.com



Each Stereodraft comes complete with:

Connected water tubing and 2 types of connectors for either a Gaymar/Stryker or Colder connector heat pump.
 Adapters: 90 degree, luer lock quick disconnect, and T.
 3 adhesive backed Velcro tubing holders.
 10' medical grade white tubing.



Model Number	Description	List
SD-M-01	Stereodraft, mouse. Sized: 3 1/2" wide x 6" long x 3/4" high	\$780.00
SD-R-01	Stereodraft, rat. Sized: 3 1/2" wide x 10" long x 3/4" high	\$805.00